

Product datasheet for RR202398L4

Spon1 (NM_172067) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spon1 (NM_172067) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Spon1
Synonyms:	Sponf
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202398).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_172067
ORF Size:	2421 bp



[View online >](#)

