

## Product datasheet for **RR202437L3V**

### Adamts1 (NM\_024400) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Adamts1 (NM_024400) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adamts1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024400
ORF Size:	2901 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202437).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024400.2</a> , <a href="#">NP_077376.2</a>
RefSeq Size:	4395 bp
RefSeq ORF:	2904 bp
Locus ID:	79252
Cytogenetics:	11q11
Gene Summary:	disintegrin and metalloprotease that may be necessary for normal kidney morphology and function [RGD, Feb 2006]



[View online »](#)