

Product datasheet for **RR202462**

Gm2a (NM_172335) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gm2a (NM_172335) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Gm2a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR202462 representing NM_172335
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGATCGCC

ATGCGTCGTGTACCGCTGCTGCTCGTGCTGGGCTTGCTGTTCTGTCGCTGGGCTTGCTGTTTCGCTGGCCCTG
TCGCCCTTCGCGCCTCATCTCGAAGCGCCCTTCCCAACTTGGTGGCTTCTCCTGGGATAACTGTGATGA
AGGAAAGGACCTGCAGTGATCAAAGCCTCACGCTCCAACCTGACCCCATTTGTCGTTCTTGAGATGTG
ATCGTCAGTGCTGAGGGCAAGACAGCATTCCCCTCACTTCTCCTCAGAAGGTGGAGCTCACCGTGGAGA
AGGAAGTGGCTGGCTTCTGGGTCAAGATCCCTTTCGCTAGAACAGCTAGGAAGCTGTACCTATGAGAATGT
CTGTGACCTGATAGACCAATACATCCCCCTGGAGAGACCTGCCAGAGCCGCTGCACACCTACGGGCTG
CCCTGCCATTGTCCCTTCAAGGAAGGCACCTACTCACTGCCTTCGAGCAACTTCACAGTGCCTGATCTGG
AGCTTCCAAGCTGGTAAGCACGGGCAACTACCGCATCCAGAGCATCTTGAGCAGCGGTGGAAGCGCCT
GGCTGCATCAAGATTGCCGCTCTCTCAAGGGCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR202462 representing NM_172335
Red=Cloning site Green=Tags(s)

MRRVPLLLVLLGLLFLVLLFAGPVAPSRLISKRPSQLGGFSWDNCDEGKDPVAVIKSLTLQDPDIVVPGDV
IVSAEGKTSIPLTSPQKVELTVEKEVAGFWVKIPVVEQLGSCYENVCDLIDQYIPPGETCEPLHTYGL
PCHCPFKEGTYSLPSSNFTVPDLELPSWLSTGNYRIQSILSSGGKRLACIKIAASLKGK

TRTRPLEQKLISEEDLAANDILDYKDDDDK

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_172335

ORF Size: 597 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_172335.2</u> , <u>NP_758838.2</u>
RefSeq Size:	2028 bp
RefSeq ORF:	600 bp
Locus ID:	282838
Cytogenetics:	10q22
MW:	21.5 kDa
Gene Summary:	mouse homolog is involved in sphingolipid metabolism [RGD, Feb 2006]