

Product datasheet for **RR202469L4V**

Dpm2 (NM_019252) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dpm2 (NM_019252) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dpm2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019252
ORF Size:	252 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202469).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019252.1 , NP_062125.1
RefSeq Size:	822 bp
RefSeq ORF:	255 bp
Locus ID:	29640
UniProt ID:	Q9Z325
Cytogenetics:	3p11
Gene Summary:	a regulatory enzyme; regulates biosynthesis of dolichol phosphate-mannose (DMP); catalytic activity is conferred by binding partner DPM1 [RGD, Feb 2006]



[View online »](#)