

## Product datasheet for **RR202478L3V**

### Colec12 (NM\_001025721) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Colec12 (NM_001025721) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Colec12
Synonyms:	MGC114569
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001025721
ORF Size:	2226 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202478).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001025721.1</a> , <a href="#">NP_001020892.1</a>
RefSeq Size:	2913 bp
RefSeq ORF:	2229 bp
Locus ID:	361289
UniProt ID:	<a href="#">Q4V885</a>
Cytogenetics:	18p13
Gene Summary:	mouse homolog is a scavenger receptor with a C-type lectin/carbohydrate recognition domain [RGD, Feb 2006]



[View online »](#)