

## Product datasheet for **RR202479L4V**

### **B2m (NM\_012512) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	B2m (NM_012512) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	B2m
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012512
ORF Size:	357 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202479).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012512.2</a> , <a href="#">NP_036644.1</a>
RefSeq Size:	1845 bp
RefSeq ORF:	360 bp
Locus ID:	24223
UniProt ID:	<a href="#">P07151</a>
Cytogenetics:	3q35
Gene Summary:	mouse homolog plays a role in the assembly and surface expression of MHC class I molecules [RGD, Feb 2006]



[View online »](#)