

Product datasheet for RR202495

Esrrb (NM_001008516) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Esrrb (NM_001008516) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Esrrb
Synonyms: Err2; Errb
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR202495 representing NM_001008516
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCGTCGGAAGACAGGCACCTGGGCTCTAGCTGCGGCTCCTTCATCAAGACGGAGCCATCTAGCCCAT
 CCTCGGGCATTGATGCCCTCAGCCACCACAGCCCAGCGGCTCGTCGGACGCCAGCGGTGGCTTTGGCAT
 GGCCCTGGGCACCCACGCCAACGGTCTGGACTCTCCGCCTATGTTCCGAGGTGCGGGGCTGGGAGGCAAC
 CCGTGTGCAAGAGCTACGAGGACTGTACTAGCGGTATCATGGAGGACTCGGCCATCAAGTGCAGTACA
 TGCTTAACGCCATCCCCAAGCGCTGTGCCTCGTGTGCGGGGACATTGCTTCTGGCTACCACTATGGAGT
 GGCTCCTGCGAGGCTTGAAGGCGTCTTCAAGAGAACCATTCAAGGAAACATCGAATACAGCTGCCCT
 GCCACCAACGAGTGTGAGATCACAAACGGAGGCGCAAGTCTGTGAGGCTGCCGGTTCATGAAATGCC
 TCAAAGTGGGGATGTGAAGGAAGGCGTGCCTTGACCGGGTGCAGGAGGCCGAGAAGTACAAGAG
 ACGGCTGGATTCGAGAACAGCCCTACCTGAGCTTACAGATTTCCCGCCTGCTAAAAAGCCATTGACT
 AAGATTGTCTCGTATCTACTGGTGGCCGAGCCGACAAGCTGTACGCTATGCCTCCCGACGATGTGCCTG
 AAGGGATATCAAGGCCCTGACCACTCTCTGTGACTTGGCAGATCGGGAGCTTGTGTTCCATTAGCTG
 GGCCAAGCACATCCAGGTTTCTCCAACCTGACACTCGGGGACCAGATGAGCCTGCTGCAGAGTGCCTGG
 ATGGAGATCCTCATCCTGGGCATCGTGTACCGCTCGCTTCCCTATGATGACAAGCTGGCATAACGGGAGG
 ACTATATCATGGATGAGGAACACTCTCGCCTGGTGGGGCTGCTGGAGCTTACCAGGCCATCTTGACGCT
 CGTACGCAGGTACAAGAAGCTCAAGGTGGAGAAGGAGTTTGTGATGCTCAAAGCCCTGGCCCTTGCC
 AACTCAGATTCAATGTACATCGAGAACCTGGAGGCTGTGAGAAGCTTCAGGACCTGCTGCATGAGGCGC
 TGCAGGACTATGAGCTGAGCCAGCGCCATGAGGAGCCACGGAGGGCGGGCAAGCTGCTGTTGACTGCTCC
 CCTGCTGCGGCAGACGGCAGCCAAAGCCGTCCAGCACTTCTACAGTGTGAAACTGCAGGGCAAGGTGCC
 ATGCACAACTCTTCTGGAGATGCTGGAGCCAAGGTG

ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR202495 representing NM_001008516
 Red=Cloning site Green=Tags(s)

MSEDRHLGSSCGSF IKTEPSSPSSGIDALSHHSPSGSSDASGGFGMALGTHANGLDSPPMFAGAGLGGN
 PCRKSYEDCTSGIMEDSAIKCEYMLNAIPKRLCLVCGDIASGYHYGVASCEACKAFFKRTIQGNIEYSCP
 ATNECEITKRRRKSCQACRFMKCLKVGMKKEGVRLDVRVGRQKYKRRLDSENSPYLSLQISPPAKKPLT
 KIVSYLLVAEPDKLYAMPPDDVPEGDIKALTTLCDLADRELVFLISWAKHIPGFSNLTLGDQMSLLQSAW
 MEILILGIVYRSLPYDDKLAYAEDYIMDEEHSRLVGLLELYRAILQLVRRYKCLKVEKEEFVMLKALALA
 NSDSMYIENLEAVQKLQDLLHEALQDYELSQLRHEEP RRAGKLLLTLPLLRQTAAKAVQHFYSVKLQGVKVP
 MHKLFLEMLEAKV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001008516

ORF Size:	1299 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001008516.2 , NP_001008516.2
RefSeq Size:	1302 bp
RefSeq ORF:	1302 bp
Locus ID:	299210
UniProt ID:	P11475
Cytogenetics:	6q31
MW:	48.3 kDa

Gene Summary:

Transcription factor that binds a canonical ESRRB recognition (ERRE) sequence 5'TCAAGGTCA-3' localized on promoter and enhancer of targets genes regulating their expression or their transcription activity (By similarity). Plays a role, in a LIF independent manner, in maintainance of self-renewal and pluripotency of embryonic and trophoblast stem cells through different signaling pathways including FGF signaling pathway and Wnt signaling pathways. Upon FGF signaling pathway activation, interacts with KDM1A by directly binding to enhancer site of ELF5 and EOMES and activating their transcription leading to self-renewal of trophoblast stem cells. Also regulates expression of multiple rod-specific genes and is required for survival of this cell type (By similarity). Plays a role as transcription factor activator of GATA6, NR0B1, POU5F1 and PERM1 (By similarity). Plays a role as transcription factor repressor of NFE2L2 transcriptional activity and ESR1 transcriptional activity (By similarity). During mitosis remains bound to a subset of interphase target genes, including pluripotency regulators, through the canonical ESRRB recognition (ERRE) sequence, leading to their transcriptional activation in early G1 phase. Can coassemble on structured DNA elements with other transcription factors like SOX2, POU5F1, KDM1A and NCOA3 to trigger ESRRB-dependent gene activation. This mechanism, in the case of SOX2 corecruitment prevents the embryonic stem cells (ESCs) to epiblast stem cells (EpiSC) transition through positive regulation of NR0B1 that inhibits the EpiSC transcriptional program. Also plays a role inner ear development by controlling expression of ion channels and transporters and in early placentation (By similarity).[UniProtKB/Swiss-Prot Function]