

Product datasheet for **RR202506L3V**

Ggcx (NM_031756) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ggcx (NM_031756) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ggcx
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031756
ORF Size:	2274 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202506).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031756.1 , NP_113944.1
RefSeq Size:	2754 bp
RefSeq ORF:	2277 bp
Locus ID:	81716
UniProt ID:	O88496
Cytogenetics:	4q31
Gene Summary:	catalyzes the posttranslational modification of glutamate to gamma-carboxyglutamate (Gla); may mediate functions of vitamin K-dependent proteins (VKDPs) [RGD, Feb 2006]



[View online »](#)