

Product datasheet for **RR202539L3V**

Amigo2 (NM_182816) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Amigo2 (NM_182816) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Amigo2
Synonyms:	Ali1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_182816
ORF Size:	1560 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202539).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_182816.2 , NP_877968.2
RefSeq Size:	3120 bp
RefSeq ORF:	1563 bp
Locus ID:	300186
UniProt ID:	Q7TNJ4
Cytogenetics:	7q36
Gene Summary:	transmembrane protein; may be involved in promoting neuronal survival [RGD, Feb 2006]



[View online »](#)