

Product datasheet for **RR202604**

Dgcr6 (NM_001107080) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dgcr6 (NM_001107080) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dgcr6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR202604 representing NM_001107080
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGACCCTATGCGGGCGCCGGGACGAGGGCGCGGACCGGGCCCGCAGCAGGAGCGGCACTACCAGC
TGCTGTCAGCGCTACAGAGCCTGGTCAAGAACTGCCAGCTCTTTTCAGCAGCGCTGTCTACACCAC
GCTCAGCGACTTAGCCCTGGCGCTGCTCGACGGAACGGTGTGGAGATCGTACAGGGCTTCTGGAGATC
CAGCACCTCACCGAGAAGAGCCTCTACAACCAACGACTGCGGCTGCAGAACGAACACAGAGTGCTCAGAC
AGACTCTGAGGCAAAAGCATCTGGAAGCCAGCAGTCCTGCCGGCCCCAACAACTTCCAGTGCTCCAGGC
AGCCCAGCAGCGTGAGCTAGAGGCTTGGAAACCCGGATCCGGGAGGAGCAGCAGGCAATGGACCGAAAAG
ATTGTCCTGGAGCTGGACCGGAAGTTGCTGACCAGCAGACACTGGAGAAGGCAGGGGTAGCTGGTT
TCTATGTGACCACTAATCCTCAGGAGCTGACGCTGCAGATGAAGTACTGGAAGTCAATCAGGAAGCTGCA
ACAGAGGGGCTGCCAAATGGAAAGCCGGCTCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR202604 representing NM_001107080
 Red=Cloning site Green=Tags(s)

MDRYAGAGDEAADRARQQERHYQLLSALQSLVKELPSSFQQLSYTTLSDLALALLDGTVFEIVQGLLEI
QHLTEKSLYNQRLRLQNEHRVLRQTLRQKHLEAQQSCRPHNLPVLQAAQQRELEALEHRIREEQQAMDRK
IVLELDRKVADQSTLEKAGVAGFYVTNPQELTLQMNLELIRKLQQRGCMGKPAL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001107080

ORF Size: 594 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001107080.2, NP_001100550.1</u>
RefSeq Size:	1029 bp
RefSeq ORF:	597 bp
Locus ID:	303794
Cytogenetics:	11q23
MW:	22.8 kDa