

## Product datasheet for RR202621L3

### Abhd10 (NM\_001123352) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Abhd10 (NM_001123352) Rat Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Abhd10
Synonyms:	RGD1308084
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202621).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_001123352
ORF Size:	891 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001123352.1</a> , <a href="#">NP_001116824.1</a>
<b>RefSeq Size:</b>	1251 bp
<b>RefSeq ORF:</b>	894 bp
<b>Locus ID:</b>	303953
<b>UniProt ID:</b>	<a href="#">Q5I0K5</a>
<b>Cytogenetics:</b>	11q21
<b>Gene Summary:</b>	Acts as an acyl-protein thioesterase that hydrolyzes fatty acids from acylated residues in proteins. Regulates the mitochondrial S-depalmitoylation of the nucleophilic active site residue of peroxiredoxin-5/PRDX5, a key antioxidant protein, therefore modulating mitochondrial antioxidant ability. Also catalyzes the deglucuronidation of mycophenolic acid acyl-glucuronide, an active metabolite of the immunosuppressant drug mycophenolate. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RR202621L3