

Product datasheet for **RR202647L3V**

Adgrl4 (NM_022294) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adgrl4 (NM_022294) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adgrl4
Synonyms:	Eltld1; Etl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022294
ORF Size:	2214 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202647).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022294.1 , NP_071630.1
RefSeq Size:	4274 bp
RefSeq ORF:	2217 bp
Locus ID:	64124
UniProt ID:	Q9ESC1
Cytogenetics:	2q45
Gene Summary:	may play a role in heart development [RGD, Feb 2006]



[View online »](#)