

Product datasheet for **RR202655L4V**

Polr1b (NM_031773) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Polr1b (NM_031773) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Polr1b
Synonyms:	RPA2; Rpo1-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031773
ORF Size:	3405 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202655).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031773.1 , NP_113961.1
RefSeq Size:	3769 bp
RefSeq ORF:	3408 bp
Locus ID:	83582
UniProt ID:	O54888
Cytogenetics:	3q36
Gene Summary:	mouse homolog is the second largest subunit of RNA polymerase I [RGD, Feb 2006]



[View online »](#)