

Product datasheet for **RR202664L4V**

PPP2CB (NM_017040) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPP2CB (NM_017040) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	PPP2CB
Synonyms:	Pp2a2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017040
ORF Size:	927 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202664).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017040.1
RefSeq Size:	1843 bp
RefSeq ORF:	930 bp
Locus ID:	24673
UniProt ID:	P62716
Cytogenetics:	16q12.3
Gene Summary:	catalytic subunit of multimeric serine/threonine phosphatase; involved in insulin signal transduction; negatively regulates the Ras/MAP kinase pathway by binding to Shc, inhibiting tyrosine phosphorylation [RGD, Feb 2006]



[View online »](#)