

Product datasheet for **RR202683L4V**

Xiap (NM_022231) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Xiap (NM_022231) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Xiap
Synonyms:	Api3; Birc4; riap3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022231
ORF Size:	1488 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202683).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022231.2 , NP_071567.1
RefSeq Size:	2469 bp
RefSeq ORF:	1491 bp
Locus ID:	63879
UniProt ID:	Q9R0I6
Cytogenetics:	Xq35
Gene Summary:	involved in directly inhibiting key apoptotic proteases, caspase 3 and 7 and is involved in maintenance of corpus luteum function and demise [RGD, Feb 2006]



[View online »](#)