

Product datasheet for **RR202688L4V**

Cyb561d2 (NM_001007753) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cyb561d2 (NM_001007753) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cyb561d2
Synonyms:	MGC95239
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001007753
ORF Size:	666 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202688).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001007753.1 , NP_001007754.1
RefSeq Size:	1500 bp
RefSeq ORF:	669 bp
Locus ID:	363137
UniProt ID:	Q641Y1
Cytogenetics:	8q32
Gene Summary:	Two-heme-containing cytochrome that catalyzes ascorbate-dependent trans-membrane ferric-chelate reduction.[UniProtKB/Swiss-Prot Function]



[View online »](#)