

Product datasheet for **RR202695L3V**

Ffar1 (NM_153304) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ffar1 (NM_153304) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ffar1
Synonyms:	Gpr40
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153304
ORF Size:	900 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202695).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153304.1
RefSeq Size:	905 bp
RefSeq ORF:	903 bp
Locus ID:	266607
UniProt ID:	Q8K3T4
Cytogenetics:	1q21
Gene Summary:	orphan G-protein coupled receptor family that is activated by medium and long chain fatty acids; may mediate physiological changes in pancreatic islets [RGD, Feb 2006]



[View online »](#)