

Product datasheet for **RR202716L4V**

Miip (NM_001017450) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Miip (NM_001017450) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Miip
Synonyms:	lip45
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001017450
ORF Size:	1158 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202716).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001017450.1 , NP_001017450.1
RefSeq Size:	1621 bp
RefSeq ORF:	1161 bp
Locus ID:	298643
UniProt ID:	Q6AYH0
Cytogenetics:	5q36
Gene Summary:	Inhibits glioma cells invasion and down-regulates adhesion- and motility-associated genes such as NFKB2 and ICAM1. Exhibits opposing effects to IGFBP2 on cell invasion (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)