

## Product datasheet for **RR202749L3V**

### Myocd (NM\_182667) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Myocd (NM_182667) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Myocd
Synonyms:	Mycd
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_182667
ORF Size:	2814 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202749).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_182667.2</a> , <a href="#">NP_872608.1</a>
RefSeq Size:	4952 bp
RefSeq ORF:	2817 bp
Locus ID:	246297
UniProt ID:	<a href="#">Q8R5I7</a>
Cytogenetics:	10q24



[View online »](#)

**Gene Summary:**

Smooth muscle cells (SM) and cardiac muscle cells-specific transcriptional factor which uses the canonical single or multiple CArG boxes DNA sequence. Acts as a cofactor of serum response factor (SRF) with the potential to modulate SRF-target genes (By similarity). Plays a crucial role in cardiogenesis and differentiation of the smooth muscle cell lineage (myogenesis) (By similarity).[UniProtKB/Swiss-Prot Function]