

Product datasheet for **RR202752L3V**

Rnaseh2b (NM_001007007) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rnaseh2b (NM_001007007) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rnaseh2b
Synonyms:	MGC94780
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001007007
ORF Size:	921 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202752).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001007007.1 , NP_001007008.1
RefSeq Size:	1450 bp
RefSeq ORF:	924 bp
Locus ID:	361056
UniProt ID:	Q5XI96
Cytogenetics:	15p12



[View online »](#)

Gene Summary:

Non catalytic subunit of RNase H2, an endonuclease that specifically degrades the RNA of RNA:DNA hybrids. Participates in DNA replication, possibly by mediating the removal of lagging-strand Okazaki fragment RNA primers during DNA replication. Mediates the excision of single ribonucleotides from DNA:RNA duplexes (By similarity).[UniProtKB/Swiss-Prot Function]