

Product datasheet for **RR202757**

Stac3 (NM_001130558) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Stac3 (NM_001130558) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Stac3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR202757 representing NM_001130558
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGACAGAAAAGGAGGTGGTGGAGTCCCCTCAGCCCCCTTCCAGGAGAGACTCCGCAAAGTGGGCTAC
AGCGACTGAAGCAATTATCAAGAAGGGCTCTCCAGAGACAACCGAGATGGAGCCTCCCCAGAGCCCCA
GGCCAATGGAGAAGCCGTGGGGGTGGGGTGGCCCCATCTACTATATCTATGAGGAAGAAGAAGAGGAG
GAGGAAGAGGAAGAACCACCCCAAGCCTCTAAGCTTGTCAATGACAAGCCCCACAAGTTCAAAGATC
ATTTCTCAAGAAGCCCAAGTTCTGTGATGTCTGTGCCGGATGATTGTGCTCAATAACAAGTTGGACT
CCGCTGTAAGAACTGCAAAACCAACATCCATGAACACTGTCACTACGTGGAGATGCAGAAATGCTTC
GGCAAGATCCCACCTGGTTTCCATCGGGCCTATAGCTCCCCACTCTACAGCAACCAGCAGTATGCTTGTG
TCTCTGCTGCCAATCGCAATGACCCTGTGTTTGAACCCTACGCGTCGGGGTGATCATGGCAAACAAGGA
ACGGAAGAAGGGCAGGCAGATAAGAAAAATCCTCTGGCAGCCATGATGGAGGAGGAGCCAGAGTCAGCC
AGGCCAGAGGAGGCAAGTCACAGGATGAAAACGATGCAGAAAAGGACAAGAAGGCTGAGAAGAAAAAC
CTGATGATAAAAAACAAGCAACCTGGTTTCCAGCAGTCTCATTACTTTGTGGCTCTCTATCGGTTCAAAGC
CCTGGAGAAGGATGATCTGGATTTCCCACCAGGGGAGAAGATCACAGTCAATTGATGACTCCAATGAGGAG
TGGTGGCGGGGAAAATCGGAGAGAAGTGGGATTCTTCCCTCAAACCTCATCATTCCGGTCCGGGCTG
GAGAACGTGTGCACCGTGAACCAAGATCCTTTGTGGGAAACCGCGAGATTGGGCAGATCACTCTCAAGAA
GGACCAGATCGTAGTGCAGAAAGGAGATGAAGCTGGTGGCTACGTCAAGGTCTACACCGGCCGAAGGTG
GGGCTGTTTCCCACCGACTTCTTGAGGAGATT

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR202757 representing NM_001130558
 Red=Cloning site Green=Tags(s)

MTEKEVVEVSPQPPFPGETPQSGLQRLKQLFKKGSPETTEMPPPEPQANGEAVGAGGGPIYYIYEEEEEE
 EEEEEPPPEPKLVNDKPHKFKDHF FKKPKFCV CARMIVLNNKFLRCKNCKTNIHEHCQSYVEMQKCF
 GKIPPGFHRAYSPLYSNQYACVSAANRNDPVFETLRVGVIMANKERKKQADKKNPLAAMEEEEPE SA
 RP EEGKSQDGNDAEKDKKAEKKT PDDKNKQPGFQQSHYFVALYRFKALEKDDLDFPPGEKITVIDDSNEE
 WWRGKI GEKVGFFPNFIIRVRAGERVHRVTRSFVGNREIGQITLKKDQIVVQKGD EAGGYVKVYTG RKY
 GLFPTDFLEEI


SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001130558

ORF Size:	1083 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001130558.1 , NP_001124030.1
RefSeq Size:	1546 bp
RefSeq ORF:	1086 bp
Locus ID:	362895
Cytogenetics:	7q22
MW:	41.1 kDa