

Product datasheet for **RR202842L3V**

Slc22a12 (NM_001034943) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc22a12 (NM_001034943) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc22a12
Synonyms:	Rst; Slc22a12; URAT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001034943
ORF Size:	1659 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202842).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001034943.1 , NP_001030115.1
RefSeq Size:	2180 bp
RefSeq ORF:	1662 bp
Locus ID:	365398
UniProt ID:	Q3ZAV1
Cytogenetics:	1q43
Gene Summary:	may act as an anion transporter; mRNA is expressed exclusively in kidney [RGD, Feb 2006]



[View online »](#)