

Product datasheet for **RR202879**

Arpc2 (NM_001106919) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Arpc2 (NM_001106919) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Arpc2
Synonyms: p34Arc
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR202879 representing NM_001106919
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATCCTGCTGGAGGTGAACAACCGCATCATCGAGGAGACGCTCGCGCTCAAGTTCGAGAACCGCGCCG
CCGAAACAACCAGAAGCTGTAGAAGTAACATTTGCAGATTTTGATGGCGTCCTTTATCATATTTCAA
TCCTAATGGAGACAAAACAAAAGTGATGGTCAGTATTTCTTTGAAATTCTACAAGAACTTCAGGCACAC
GGAGCTGATGAGTTACTAAAGAGAGTATACGGAAGTTTCTTGGTGAACCCAGAATCAGGATACAATGTGT
CTTTGCTATATGACCTGGAAAATCTACCTGCATCCAAGGATTCTATTGTGCATCAGGCTGGCATGTTGAA
GCGAAATTGTTTTGCCTCTGTGTTTGGAAAATACTTCCAATTCCAAGAAGAGGGCAAAGAAGGAGAGAAC
AGAGCAGTTATCCATTATAGAGATGATGAGACCATGATGTGGAATCTAAAAAGACAGAGTCACAGTAG
TCTTCAGCACAGTTTTTAAGGATGACGACGATGTGGTCATTGGAAAGGTGTTTCATGCAGGAGTTCAAAGA
AGGACGAAGAGCCAGCCACACAGCCCCCAGGTCCTTTTTAGCCACAGGGAACCTCCTCTAGAAGTAAA
GACACAGATGCTGCCGTGGTGACAACATTGGCTACATCACCTTTGTGCTGTTCCCTCGCCACACCAATG
CCACTGCTCGAGACAACCCATCAACTGATCCACACGTTCCGGGACTATCTGCACTACCCATTAAGTG
CTCTAAGGCCTATATTCATACACGAATGCGAGCAAAAACATCCGACTTCCTTAAGGTGCTCAACCCGTGCA
CGCCAGATGCCGAGAAAAAAGAAATGAAAACAATCACGGGAAGACTTTTTTCATCCCGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR202879 representing NM_001106919
 Red=Cloning site Green=Tags(s)

MILLEVNNRIIEETLALKFENAAAGNKPEAVEVTFADFDGVLHYHISNPNGDKTKVMVSI SLKFYKELQAH
 GADELLKR VYGSFLVNPESGYNSLLYDLENLPASKDSIVHQAGMLKRNCFA SVFEKYFQFQEEGKEGEN
 RAVIHYRDEETMYVESKKDRVT VVFSTVFKDDDDV VIGKVF MQEFKEGRRASHTAPQVLF SHREPPLELK
 DTDAAVGDNIGYITFVLFPRHTNATARDNTINLIHTFRDYLHYHIKCSKAYIHTRMRAKTSDFLKVLNRA
 RPDAEKEMKTI TGTGTFSSR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106919

ORF Size: 900 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001106919.1 , NP_001100389.1
RefSeq Size:	1428 bp
RefSeq ORF:	903 bp
Locus ID:	301511
Cytogenetics:	9q33
MW:	34.3 kDa