

## Product datasheet for **RR202903**

### Eno2 (NM\_139325) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Eno2 (NM_139325) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Eno2
Synonyms:	NSE; RNEN3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR202903 representing NM_139325 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTCTATACAGAAGATCTGGGCCGAGAGATCTTGGACTCCCGTGGGAATCCCACCGTGGAGGTGGATC  
TCCATACTGCCAAAGTCTTTCCGGGCTGCAGTCCCCAGTGGGGCTCCACTGGCATCTATGAGGCCCT  
GGAGCTAAGGGATGGGGACAAACAGCGTTACTTAGGCAAAGGTGTCTGAAGGCTGTGGACCACATCAAC  
AGCACCATCGCACCGGCCCTCATCAGCTCAGGCCCTCTCTGTGGTGGAGCAGGAGAAGCTGGACAACCTGA  
TGCTGGAGTTGGATGGGACTGAGAACAATCCAAGTTTGGGGCCAATGCCATCCTGGGTGTGTCCCTGGC  
CGTGTGCAAGGCTGGGGCAGCCGAGAAGGACTTCCCCCTATCGCCACATTGCTCAACTGGCTGGGAAC  
TCCGACCTCATCTGCCGTGCCGGCCTTAATGTGATCAACGGTGGCTCTCATGCTGGGAACAAGTTGG  
CCATGCAGGAGTTCATGATCCTCCAGTGGGTGCTGAGAGCTTTCGGGATGCCATGCGACTTGGGGCCGA  
GGTGTACCACACACTCAAGGGGGTCAAGGACAAGTACGGCAAGGATGCCACTAATGTGGGGATGAA  
GGCGGCTTTGCCCCAAATCCTGGAGAACAGCGAAGCTTGGAGCTGGTGAAGGAAGCCATTGACAAGG  
CTGGCTACACGGAAAAGATGGTATTGGTATGGATGTGGCTGCCCTGAGTTTTACC GCGATGGCAAATA  
CGACTTGGATTTCAAGTCTCCCGCTGACCTTCCCGATGCATCACTGGGACCAGCTTGGGGCACTCTAC  
CAGGACTTTGTCCGAACTATCCTGTGGTCTCCATTGAAGACCCATTCCGACCAGGATGACTGGGCAGCTT  
GGTCCAAGTTCACAGCCAATGTCGGCATCCAGATAGTGGGTGATGACCTGACGGTGACCAACCCCAAGCG  
CATCGAGCGGGCAGTGGAGGAGAAGGCCTGCAACTGTTTGTCTCAAGGTCAACCAGATCGGCTCAGTC  
ACAGAAGCCATCCAAGCGTGCAAGCTGGCCAGGAGAACGGCTGGGGGTTATGGTGAGTCATCGCTCTG  
GAGAAACCGAGGACACGTTTCATTGCAGACCTCGTAGTGGGACTGTGTACAGGTGAGTCAAGACTGGTGC  
CCCATGCAGATCTGAACGTCTGGCGAAGTACAACCAGCTCATGAGGATTGAAGAGGAGTTGGGGAGGAG  
GCTCGCTTCGCGGACACAACCTCCGGAATCCCAGTGTGCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



**Protein Sequence:** >RR202903 representing NM\_139325  
Red=Cloning site Green=Tags(s)


MSIQKIWAREILDSRGNPTVEVDLHTAKGLFRAAVPSGASTGIYEALERDGDQRYLKGKVLKAVDHIN  
 STIAPALISSGLSVVEQEKLNDLMELDGTENKSKFGANAAILGVSLAVCKAGAAEKDPL YRHIAQLAGN  
 SDLILPVPAFNVINGGSHAGNKLAMQEFMILPVGAESEFRDAMRLGAEVYHTLKGVIKDKYGDATNVGDE  
 GGFAPNILENSEALELVKEAIDKAGYTEKMGVIGMDVAASEFYRDGKYDLDFKSPADPSRCITGDQLGALY  
 QDFVRNYPVVSIEDPFDQDDWAASKFTANVGIQIVGDDLTVTNPKRIERA VEEKACNCLLLKVNQIGSV  
 TEAIQACKLAQENGWGMVSHRSGE TEDTFIADLVVGLCTGQIKTGAPCRSERLAKYNQLMRIE EELGEE  
 ARFAGHNFRNPSVL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:**

NM\_139325

<b>ORF Size:</b>	1302 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_139325.3</a> , <a href="#">NP_647541.1</a>
<b>RefSeq Size:</b>	2222 bp
<b>RefSeq ORF:</b>	1305 bp
<b>Locus ID:</b>	24334
<b>UniProt ID:</b>	<a href="#">P07323</a>
<b>Cytogenetics:</b>	4q42
<b>MW:</b>	47.1 kDa
<b>Gene Summary:</b>	catalyzes the interconversion of 2-phosphoglycerate to phosphoenolpyruvate; may play a role in glycolysis [RGD, Feb 2006]