

Product datasheet for **RR202986L3V**

Cpt1c (NM_001034925) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cpt1c (NM_001034925) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cpt1c
Synonyms:	CPT1-B; CPTI-B; CPT IC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001034925
ORF Size:	2403 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202986).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001034925.2 , NP_001030097.2
RefSeq Size:	2747 bp
RefSeq ORF:	2406 bp
Locus ID:	308579
UniProt ID:	F1LN46
Cytogenetics:	1q22
Gene Summary:	May play a role in lipid metabolic process.[UniProtKB/Swiss-Prot Function]



[View online »](#)