

Product datasheet for **RR202993**

Tmem184a (NM_001025413) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem184a (NM_001025413) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tmem184a
Synonyms:	RGD1306702
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR202993 representing NM_001025413 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGGAATGCATCAGGGTTTCTGAAGACAGCTGGAGCCCCCTTGTGTCAGCGACCTGGCTTCGCCCA
GCCCCCTCCGGCAATGCCATGGTGGCAGCCGGCCACAGATGGAGCGCGTGGACAACGGCTCGCAGGG
GGCCCCCAGCTCTTCTCACCAGCGCATTGGCCCGAGGGCTTTCAGGCGTATTTGTATGGACTGCTCTG
CTGCTAACCGGCCACCAGATCTACTCCCACCTGGGTTCTTACACCGTCCCAGGAGAGCAGCGCTTGTCA
TCCGCCTCCTATTGTCGCCATCTACGCCTTCGACTCATGGCTCAGCCTCCTCCTCGGGGGCCA
TCCTTACTACGTCTACTTCGACTCTGTGCGAGACTGCTACGAAGCATTGTGTCATCTACAGTTTCTGACC
TTGTGCTTCCAATATCTGGGGGGCAGAGCGCCATCATGGCTGAGATCAGAGGCAAACCTATCAGGTCCA
GCTGCTTCTATGGGACCTGCTGCCTCCGTGGTATGTCCTACTCCATCACTTTCTTACGCTTCTGCAACA
GGCCACACTCCAGTTCTGCATTGTGAAACCCGTTATGGCGCTGATCACCATCATCCTCCAGGCTTTTGAC
AAATACCAGATGGCGACTTCAACATCCACAGCGGCTACCTGTACGTGACCCTCGTGTACAACGCCTCTG
TTAGCCTGGCTCTCAGCTCTGTTCTTTCTACTTCGCTACCAGGGACCTCCTGAGACATTGAGCC
AGTACTCAAGTTCTCACCATCAAAGCCATCATCTTCTCCTTCTGGCAAGGAATGCTGTTGGCCATC
TTGGAGAGGTGTGGGTCATCCCTGAGGTCCAGGCCGTGGACGGCACCAGGGTTGGGGCTGTTACCCTAG
CCGCTGGCTACCAGAATTTCTCATCTGTATTGAGATGCTGTTTGCCTCCCTTGGCTGCGTTACGCCTT
CCCCAGCCAGGTGATTAGAAAAGAAGAACTCACCAGCCCCACCAGCTCCTATGCAAAGCATATCCAGT
GGCTTAAAGGAGACCATCAGCCACAGGACATTGTCCAGGATGCCATTACAACCTTCTCACCAGCCTACC
AGCAGTATACACAGCAGTCTACACATGAAGCCCTGGGCTGGCCAGGGCGGGCACCATCGCCAGCAC
CCACCTGGCCAGCCAGTGGCTCTGGAGGTGGTAAGAAGAGTCGCAACATAGAGAAACGCATGCTGATC
CCCTCAGAGGACCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR202993 representing NM_001025413
 Red=Cloning site Green=Tags(s)


MRNASGFLKTAGAPLVSATWLPSPPPAMPVMAAGPQMERVDNGSQGAPQLFLTSALARGVSGVFVWTA
 LLTGHQIYSHLRSYTVPREQRFVIRLLFIVPIYAFDSWLSLLLLGGHPYVYFDSVRDCYEFVIYSFLT
 LCFQYLLGGESAIMAEIRGKPIRSCFYGTCCLRGMSYSITFLRFCKQATLQFCIVKPMALITIIILQAFD
 KYHDGDFNIHSGYLYVTLVYNASVSLALYALFLFYFATRDLRPFEPVLKFLTIIKAIIFLSFWQGMLLAI
 LERCGVIPEVQAVDGRVAGATLAAGYQNFLICIEMLFASLALRYAFPSQVYSEKKNPAPPAPMQSISS
 GLKETISPDIVQDAIHNFSPAYQQYTQSTHEAPGPGQGGHPSPSTHPGPASGSGGGKSRNIEKRMLI
 PSEDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001025413

ORF Size:	1275 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001025413.1 , NP_001020584.1
RefSeq Size:	1855 bp
RefSeq ORF:	1278 bp
Locus ID:	304325
UniProt ID:	Q4QQS1
Cytogenetics:	12q11
MW:	46.9 kDa
Gene Summary:	Acts as a heparin receptor in vascular cells (PubMed:26769966). May be involved in vesicle transport in exocrine cells and Sertoli cells (By similarity).[UniProtKB/Swiss-Prot Function]