

Product datasheet for **RR203009L4V**

Snx24 (NM_001008364) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Snx24 (NM_001008364) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Snx24
Synonyms:	MGC94664
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001008364
ORF Size:	507 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203009).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008364.1 , NP_001008365.1
RefSeq Size:	2216 bp
RefSeq ORF:	510 bp
Locus ID:	361328
UniProt ID:	Q5U2S5
Cytogenetics:	18q11
Gene Summary:	May be involved in several stages of intracellular trafficking.[UniProtKB/Swiss-Prot Function]



[View online »](#)