

Product datasheet for **RR203020L4V**

Ece2 (NM_001002815) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ece2 (NM_001002815) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ece2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001002815
ORF Size:	2295 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203020).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001002815.2 , NP_001002815.2
RefSeq Size:	3010 bp
RefSeq ORF:	2298 bp
Locus ID:	408243
Cytogenetics:	11q23
Gene Summary:	metalloprotease type II enzyme; catalyzes functionally active forms of the endothelin vasoactive peptide family [RGD, Feb 2006]



[View online »](#)