

Product datasheet for **RR203022**

Rfc4 (NM_001105869) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rfc4 (NM_001105869) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rfc4
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR203022 representing NM_001105869 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAAGCATTCTGAAAGGCACATCTGTCTAGTACCAAACCACAATTGACCAAGGAACGGGGCATATCTT
CCTCTGCAGGAAGCAGTGGAGAAACCAAGAAAGTCAAACCAGTTCCTTGGGTGAAAAATACCGTCCGAA
GTGCGTGGATGAAGTTGCTTTCCAGGAAGAAGTGTTCAGTGCCTGAAAAAGTCTTTAGAAGGAGCTGAT
CTTCCTAATCTCTTGTCTATGGCCACCTGGAAGTGGCAAACATCCACCATTTTGGCAGCAGCTAGAG
AGCTCTTTGGGCCTGAACTCTTTCGATTAAGAGTCTTGAAGTAAATGCATCTGATGAACGTGGAATACA
AGTAGTCCGAGAGAAAGTAAAAACTTTGCTCAGTTAACTGTGTCTGGAAGTCGTTTCCAGTGGAAACCA
TGTCTCCCTTTAAGATTGTAATTCTGGATGAAGCAGATTCTATGACTTCAGCTGCTCAGGCAGCTTTAA
GACGAACCATGGAAAAAGAGTCTAAAACAACCAGATTCTGCCTCATCTGTAACATGTCAGCCGAATAAT
TGAACCCCTCACTTCTAGATGTTCAAAATTCGCTTCAAGCCTCTGTGAGATAAAATTCAACAGAAGCGA
TTATTGGATATTGCTGAGAAGGAAAATGTCAAAATTTGGCGATGAGGAAAAGCTTATCTTGTGTAATAT
CAGAAGGAGATTTACGGAAAGCAATTACATTTCTCAAAGTGTACTCGACTAACAGGTGGAAGGAAAAT
CTCGGAAGATGTGATCACAGACATTGCTGGGGTAAATACCAGCTGCAACCATTTGAGGGAATAGTCACTGCA
TGTACAGCGGCTCTTTTGACAACTAGAAGCCGTGCTGAAGAACCCTAATAGATGAAGGGCATGCAGCCA
CTCAGCTTGTCAACCACTTCATGATTCATCATAGAAGATGAAAATCTGTCTGATAAACAGAAGTCCAT
CATTACAGAAAACTTGCTGAGGTGGATAAGTGTTTAGCAGATGGTGCAGATGAGCACCTGCAGCTGATG
AGCCTCTGCGCAACTGTGATGCAGCAACTGACCCAGAAGTGT

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR203022 representing NM_001105869
 Red=Cloning site Green=Tags(s)

MQAFLKGTSVSTKPQLTKERGISSAGSSGETKKVKPVPVWEKYRPKCVDEVAFQEEVAVLKKSLLEGAD
 LPNLLFYGPPGTGKTSTILAAARELFGPELFLRVLELNASDERGIQVREKVKNF AQLTVSGSRSDGKP
 CPPFKIVILDEADSMTSAAQAALRRTMEKESKTRFCLICNYVSRIEPLTSRCSKFRFKPLSDKIQQKR
 LLDIAEKENVKIGDEEIAYLVRISSEGLRKAITFLQSATRLTGGKEISEDVITDIAGVIPAAATIEGIVTA
 CHSGSFDKLEAVLKNLIDEGHAATQLVNQLHDSIIEDENLSDKQKSIITEKLAEVDKCLADGADEHLQLM
 SLCATVMQQLTQNC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001105869

ORF Size:	1092 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001105869.2 , NP_001099339.1
RefSeq Size:	1265 bp
RefSeq ORF:	1095 bp
Locus ID:	288003
Cytogenetics:	11q23
MW:	40 kDa