

Product datasheet for **RR203040L3V**

Sub1 (NM_001009618) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sub1 (NM_001009618) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sub1
Synonyms:	Rpo2tc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001009618
ORF Size:	381 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203040).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009618.1 , NP_001009618.1
RefSeq Size:	782 bp
RefSeq ORF:	384 bp
Locus ID:	192269
UniProt ID:	Q63396
Cytogenetics:	2q16
Gene Summary:	expressed in 12 day embryo and in a chemically induced pancreatic B-cell tumor; not expressed in normal differentiated cells [RGD, Feb 2006]



[View online »](#)