

## Product datasheet for **RR203051**

### Arhg dib (NM\_001009600) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Arhg dib (NM\_001009600) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Arhg dib  
**Synonyms:** MGC108926  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR203051 representing NM\_001009600  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGACGGAGAAGGATGAGCAGCCACAGCTGGAGGAGGGAGAGGATGAGCTGGACAGCAAGCTCAACTACA  
AGCCGCCCCCTCAGAAGTCCCTGAAGGAGCTGCAGGAGATGGACAAAGATGACGAGAGTCTAACCAAGTA  
CAAGAAAACACTGCTGGGAGATGTCCCGTGGTAGCAGACCCAACAGTCCCAACGTGACTGTGACCCGG  
CTTAGCCTTGTGGTGTGACAGTGACACAGGACCAATCACCATGGACCTTACTGGTGTCTCAAGGCGCTCA  
AAAAGGATACATTTGTGCTAAAGGAAGGTATCGAATACAGGGTCAAGATTAACCTCAAAGTGAATAAGGA  
TATTGTGTCTGGCCTGAAGTATGTTTACGACAGACGTACCGGAATGGGATGAAAGTGGATAAAGCTACGTTC  
ATGGTTGGCAGCTACGGGCTCGACCAGAGGAGTACGAGTTTCTCACTCCAGTAGAGGAGGCTCCCAAGG  
GCATGTTGGCCCGAGGCACTTACCACAACAAGTCTTCTTACCAGATGATGACAAGCAGGACCACCTCAC  
CTGGGAATGGAACCTGGCCATTAAGAAGGATTGGACAGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR203051 representing NM\_001009600  
Red=Cloning site Green=Tags(s)

MTEKDEQPQLEEGEDELDSKLNYPKPPQKSLKELQEMDKDDESLTKYKKTLLGDVPPVADPTVPNVTVTR  
LSLVCD SAPGPI TMDL TGD LKALKKDTFVLKEGIEYRVKINFKV NKDIVSGLKYVQQT YRNGMKV DKATF  
MVGSYGPRPEEYEF LTPVEEAPKGM LARGTYHNKSF FTDDDKQDHL TWENLAIKKDWT D


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_001009600

**ORF Size:** 600 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** NM\_001009600.1, NP\_001009600.1

**RefSeq Size:** 1157 bp

**RefSeq ORF:** 603 bp

**Locus ID:** 362456

**Cytogenetics:** 4q43

**MW:** 22.9 kDa