

Product datasheet for **RR203068L4V**

Tbx3 (NM_181638) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tbx3 (NM_181638) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tbx3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_181638
ORF Size:	2229 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203068).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181638.1 , NP_853669.1
RefSeq Size:	2232 bp
RefSeq ORF:	2232 bp
Locus ID:	353305
UniProt ID:	Q7TST9
Cytogenetics:	12q16
Gene Summary:	human homolog acts as a transcriptional repressor and plays a role in the development of several organ systems [RGD, Feb 2006]



[View online »](#)