

## Product datasheet for **RR203071**

### Pofut2 (NM\_001107621) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Pofut2 (NM_001107621) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pofut2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR203071 representing NM_001107621 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGGCTCAGCGTCGTCTGTCTGCTGCTGGCAGCAGCATCCTGGCGCCGGTCTCGGCTTCAGGAG  
AAGAATTCTGGCCTGGCCAGTCGGCGCCGACATCCTGTCCGGGGCGGCATCCCGCAGACGGTACCTTCT  
GTATGATGTCAACCCCCGGAAGGCTTTAATCTCCGCAGGGATGTCTATATCCGGTGCCATCCTTGCTA  
AAGACGCTGCTGAAGACAGAGGAGTGGGTGCTTGTCTGCCCCATGGGGCCGCCTTACCAGTGGCAAA  
GCCCTGATATCCATCAGGTCCGGATTCCCTGGTCTGAGTTTTTTGACCTCCCGAGTCTCAATAAAAAACAT  
CCCCGTTATAGAGTACGAGCAGTTCATTGCAGAGTCCGGCGGGCCCTTTATTGACCAGGTGTACGTCCTG  
CAAGGCTACGCCGAGGGCTGGAAGGAGGGCACCTGGGAGGAGAAGGTGGACGAGCGGCCCTGCATAGACC  
CTCTGCTGTACTCACAGGACAAGCATGAGTATTTAGAGGCTGGTTTTGGGGCTACGAAGAAACACGGGG  
TCTGAATGTCTCCTGCCTGTCCGTCCAGGGTCCGGCTTCTGTATCGCGCCTGTACTCTAAAAAACACA  
TCTGCTCGGTCTGTGATGCTGGACCGAGCTGAAAACCTGCTTCATGACCACTATGGAGGCAGAGAGTACT  
GGGACACCCGGCGCAGCATGGTGTGGCAAAGCACTTGGCGGCCGTGGGAGACGAATTCAGGAGCCAGCA  
CCTCAACTCCACGGATGCCGCCGATAAGATGGTCTGGAGGAAGACTGGACCAGGATGAAGGTCAAAGT  
GGCTCAGCGCTAGCGGCCCTACCTTGGCGTCCACCTGAGAAGGAAGGACTTCACTGGGGCCACAGGG  
AGGATGTGCCAGCCTAGAAGGTGCCGTGAAGAAGATCCGCAGCCTCATGAAGACTCACCAGCTGGACAA  
GGTGTTCGTGGCCACAGATGCCATCAGGAAGGAACAGGAAGAAGTGAAGGAGCTGCTGCCGAAATGGTG  
AGGTTTGAAGCCAGTGGGAGGAGCTGGAAGTGTACAAGGACGGAGGCGTCGCCATCATCGACCAGTGGA  
TCTGCGCTCATGCCAGGTTTTTTATTGGTACCTCCGTTTCCACGTTTTCCTTTCGGATTCATGAAGAGAG  
AGAGATCTGGGGTTGGACCCCAAGACGACATAACAACCGTTCGTGGAGACCAGGAGAAAGCGTGCCGAG  
CAGCCCACGCACTGGAAGATTGCGTAC

**ACGCGT**ACGCGGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR203071 representing NM\_001107621  
 Red=Cloning site Green=Tags(s)


MAALSVVCLLLAAASWRPVSASGEEFWPGQSAADILSGAASRRRYLLYDVPPEGFNLRDQVYIRVASLL  
 KTLKTEEWLVLPPWGRLYHWQSPDIHQVRIPWSEFFDLPSLNKNIPVIEYEQFIAESGGPFIDQYVVL  
 QGYAEGWKEGTWEEKVDERPCIDPLLYSQDKHEYFRGWFWGYEETRGLNVSVLSVQGSASVIAPVLLKNT  
 SARSVMLDRAENLLHDHYGGREYWDTRRSMVFAKHLRAVGDEFRSQHLNSTDAADKMVLEEDWTRMKVKL  
 GSALGGPYLGVHLRRKDFIWGHREDVPSLEGAVKKIRSLMKTHQLDKVVFATDAIRKEQEELRKLPEMV  
 RFEPTWEELELYKDGVAIIDQWICAHARFFIGTSVSTFSFRIHEEREILGLDPKTTYNRFCDGQEKACE  
 QPTHWKIAY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



<b>ACCN:</b>	NM_001107621
<b>ORF Size:</b>	1287 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001107621.2</a> , <a href="#">NP_001101091.2</a>
<b>RefSeq Size:</b>	2259 bp
<b>RefSeq ORF:</b>	1290 bp
<b>Locus ID:</b>	309686
<b>Cytogenetics:</b>	20p12
<b>MW:</b>	49.5 kDa