

Product datasheet for **RR203075**

Socs4 (NM_001107256) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Socs4 (NM_001107256) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Socs4
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR203075 representing NM_001107256 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGAAAACAATAGTAAAAATGTAGATGTGCGGCCTAAAACAAGTCGGAGTAGAAGTGCTGACAGGA
AGGACGGTTATGTGTGGAGTGAAAGAAGTTATCTTGGTCCAGAAAGAGTGAGAGCTGTTCTGAGTCTGA
GGCCGTAGGTACTCTTGAGAGCGTTGAAATCCCTCTGAGGAGCCAGGAAAGGGAACCTAGCTGTTCTGTC
ATTGAGTTGGACTTGGATCATTCCCTGTGGGCATAGATTTTAGGCCGATCCCTTAAACAGAAACTGCAAG
ATGCCGTGGGGCAGTGTTTTCTATAAGGAACTGTGGTGCTCGCCACTCTCCAGGGCTTCCATCTAAAAG
AACGATTCACATCAGTGAGCTCATGCTAGACAAGTGCCCTTTCCACCTCCATCAGACTTAGCCTTTAGG
TGGCATTATTAACGGCACACTGTTGCTCTGAGTTCAGTTCAGATGAATGGGTGAGCACAGACCTGT
GTGAGAGGAGGCCGAGAGATGCTCAGCTGAAACGAAGAAGCACAGAAGACGACATACCCTGTTTCTCACA
TACCAGTGCCAGCCGTGTGCATAACTACCAACAGTGCTTCATGTGCAGGTGGTCATGTGACCGGTTCT
CTGATGAACTTGGTCACCGACAATAGCATGGAAGACGGTGATATGGACTCAGAGGATGAAATTATAACAC
TGTGCACAAGCTCCAGGAAAAGGAACAAGCCAGATGGGAAACGGATGACGGCATCCTGCAATTGGAGGC
ACCTCCTAAGTTCCACACCCAGATTGATTACGTGCACTGCCTTGCCAGACCTCCTCCAGATCAGTAAC
AACCCGTGCTACTGGGGGTTATGGACAATACGCAGCCGAGGCTCTGCTAGAAGGAAAGCCAGAGGGCA
CCTTTTACTTCGAGATTCAGCCAGGAAGATTATTTATCTCTGTTAGTTTTAGGCGCTACAGTCGTTT
TCTTACGCTAGAATTGAACAGTGGAATCATAACTTCAGCTTTGACGCCACGATCCCTGTGTCTCCAT
TCTCCGATATTACTGGGCTCCTAGAGCATTATAAGGACCAAGTGCCATGTTCTTTGAGCCACTTT
TATCCACTCCGTTAATCCGGACTTTCCCTTTTCTTGCAGCACATTTGCAGGACAGTTATTTGTAAGT
TACAACGTACGATGGCATCGATGCCCTTCAATTCCTTCTCTATGAAATTGTACCTGAAGGAATACCAT
TACAAATCAAAGTTAGGCTACTCAGGATTGACGTGCCAGAGCCACAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR203075 representing NM_001107256
 Red=Cloning site Green=Tags(s)


MAENNSKNVDVRPKTSRSRSADRKDGIVWSGKLSWSRKSESCSESEAVGTLESVEIPLRSQERELSCSS
 IELDLDHSCGHRFLGRSLKQKLQDAVGQCFPIRNCGARHSPGLPSKRTIHISELMLDKCPFPSPDLAFR
 WHFIKRHTVALSSSSDEWVSTDLRERRRDAQLKRRSTEDDIPCFSHTSVQPCVITTNASACAGGHVTGS
 LMNLVTDNSMEDGMDSEDEIITLCTSSRKRKPRWETDDGILQLEAPPKFHTQIDYVHCLVPDLLQISN
 NPCYWGVMKYAAEALLEGKPEGTFLLRDSAQEDYLFVSVFRRYSRSLHARIEQWNHNFSDAHDPCVFH
 SPDITGLLEHYKDPACMFFPELLSTPLIRTFPFLQHCRTVICNCTTYDGIDALPIPSMKLYLKEYH
 YKSKVRLLRIDVPEPQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001107256
ORF Size:	1308 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107256.2 , NP_001100726.1
RefSeq Size:	2037 bp
RefSeq ORF:	1311 bp
Locus ID:	305828
Cytogenetics:	15p14
MW:	49.6 kDa