

Product datasheet for **RR203080L4V**

Taar1 (NM_134328) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Taar1 (NM_134328) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Taar1
Synonyms:	taR-1; Tar1; Trar1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_134328
ORF Size:	996 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203080).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134328.1 , NP_599155.1
RefSeq Size:	999 bp
RefSeq ORF:	999 bp
Locus ID:	113914
UniProt ID:	Q923Y9
Cytogenetics:	1p12
Gene Summary:	G-protein coupled receptor that induces cAMP production in response to para-tyramine and other trace amines; also acts as a receptor for catecholamine metabolites and amphetamines [RGD, Feb 2006]



[View online »](#)