

Product datasheet for **RR203083**

Ube2j1 (NM_001106642) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ube2j1 (NM_001106642) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ube2j1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR203083 representing NM_001106642
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGACCCGCTATAACCTGAAGAGCCCGCTGTAAACGATTGATGAAAGAAGCCGAGAATTGAAAG
ATCCAACAGATCATTACCACGCTCAGCCTCTGGAGGATAACCTTTTTGAATGGCACTTCACAGTTAGAGG
GCCCCAGATCCGATTTTATGGAGGAGTTTATCATGGACGCATAGTACTGCCACCAGAGTACCCCATG
AAACCACCAGCATTATTCTGCTGACGGCTAATGGACGATTTGAAGTGGGCAAGAAGATCTGTTTAAAGCA
TCTCAGGACATCATCCTGAGACCTGGCAGCCTTTCATGGAGCATAAGGACAGCCCTGCTGGCCATCATTGG
GTTTATGCCAACTAAAGGAGAGGGAGCCATAGGTTCTCTAGACTACACGCCCGAGGAGAGGAGAGCACTT
GCCAAGAAATCACAAGACTTCTGTTGTGAGGGCTGTGGCTCTGCCATGAAGGATGCTCTGTTGCCTTTAA
AATCTGGAATCGATTCAAGCCAGGCTGACCAAGAAGCCAAAGAACTGGCTAGACAAATTAGTTTTAAGGC
CGAAGTCAATTCATCTGGAAAGACTATCGCTGAATCAGACTTAAACCAGTCTTCTCACTAAATGATTCA
CAAGACGATTACCTACAACATCCAGGGTGTGCAGCAAGCACATCGTATGGCGCCAGAATCCATCAGG
AGCACCTTGCCACAACCTACCCAGCCAGCTCCTAAGAATACCTCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR203083 representing NM_001106642
Red=Cloning site Green=Tags(s)


METRYNLKSPAVKRLMKEAAELKDPTDHYHAQPLEDNLFEWHFTVRGPPDSDFDGGVYHGRIVLPPEYPM
KPPSIILLTANGRFVGGKICLSISGHPETWQPSWSIRTALLAIIGFMPTKGEAIGSLDYTPERRAL
AKKSQDFCCEGCGSAMKDVLLPLKSGIDSSQADQEAKEALARQISFKAENVSSGKTI AESDLNQSFSLNDS
QDDYLQHSRVLQQAHRMAPRIHQEHPCHNLPSQLLRIPP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV




Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106642

ORF Size: 747 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106642.1, NP_001100112.1</u>
RefSeq Size:	1306 bp
RefSeq ORF:	750 bp
Locus ID:	297961
Cytogenetics:	5q21
MW:	27.9 kDa