

## Product datasheet for **RR203090**

### Mpzl1 (NM\_001007728) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mpzl1 (NM\_001007728) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Mpzl1  
**Synonyms:** MGC93692  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR203090 representing NM\_001007728  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCAGAGGCCGTCGGAGCCGTGACATTGATTGCGGCCCGGCCCGCCGGCGTGGCTGTGGTCGGCGC  
TAGCCCGGATGCTCGGGCTGTTGACAGCTAGAATATCAGCCTTGGAGGTACATACTCTAAAGAAATCTT  
TGTGGTGAATGGGACACAAGGGAAGCTCACGTGCACGTTTGACTCTCCTAATACTGGGTGGTTGACC  
ACCGTCTCCTGGAGCTTCCAGCCAGAAGGCACCGACAGCGCCGTGTCGTTTTTCCACTACTCACAAGGAC  
AAGTGTACATCGGGATTATCCACCGTTTAAAGACCGAGTCACCTGGGCTGGGGACCTTGACAAGAAAGA  
TGCATCGATCAATATAGAAAACATTCAGGCTGTTCAATGGCACCTATATCTGCGATGTCAAAAACCTT  
CCTGACATTGTTGTGCGGCTGGGCAGATTAGGCTCCAATGGTGGAAATAGACAACCTGCTGGTGTTC  
TGGTTTGGGTGGTGGGCACTGTCACTGCTGTGGTCTTGGGCTCACTCTACTCATCAGCTTGGTCTT  
GGTCGCTCTACAGACGGAAGCATTGAAGCGGGATTATACCGGCTGCAGTACCTCAGAGCGTTTGTCA  
CCAGTTAAGCAGGCTCCACGGAAGTGTCCCTCCGACACAGAGGGTCTAGTAAAGAGTCCGCTTCCGCCG  
GATCTCACCAGGCCCAGTCATATACGCACAGTTAGACCACTTGCGGGACATCACAGTGGAAGATTAA  
TAAGTCAGAGTCTGTTGTATGCGGACATCCGAAAGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR203090 representing NM\_001007728  
 Red=Cloning site Green=Tags(s)

MAEAVGAVTLIAAPARRRWLSALAAMLGLLTARISALEVHTPK EIFVVTNGTQGKLTCTFDSPNTTGWLT  
 TVSWSFQPEGTD SAVSFFHYSQGVYIGDYPPFKDRVTWAGDLKKDASINIENIQAVHNGTYICDVKNP  
 PDIVVVRPGQIRLHVVEIDNLLVFLVWVVVGTAVVLGLTLLISLVLVLYRRKHSKRDYTGCTSERLS  
 PVKQAPRKCPSDTEGLVKSPPSAGSHQGPVIYAQLDHSGGHHSKGINKSESVVYADIRKD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001007728  
 ORF Size: 810 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001007728.1</a> , <a href="#">NP_001007729.1</a>
<b>RefSeq Size:</b>	2291 bp
<b>RefSeq ORF:</b>	813 bp
<b>Locus ID:</b>	360871
<b>UniProt ID:</b>	<a href="#">Q6AYT8</a>
<b>Cytogenetics:</b>	13q23
<b>MW:</b>	29.5 kDa
<b>Gene Summary:</b>	Cell surface receptor, which is involved in signal transduction processes. Recruits PTPN11/SHP-2 to the cell membrane and is a putative substrate of PTPN11/SHP-2. Is a major receptor for concanavalin-A (ConA) and is involved in cellular signaling induced by ConA, which probably includes Src family tyrosine-protein kinases. May be involved in regulation of integrin-mediated cell motility (By similarity).[UniProtKB/Swiss-Prot Function]