

Product datasheet for **RR203117L4V**

Rem2 (NM_022685) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rem2 (NM_022685) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rem2
Synonyms:	MGC156645
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022685
ORF Size:	1023 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203117).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022685.2 , NP_073176.2
RefSeq Size:	1835 bp
RefSeq ORF:	1026 bp
Locus ID:	64626
UniProt ID:	Q9WTY2
Cytogenetics:	15p13
Gene Summary:	a Ras-related small GTP-binding protein that is found in high levels in neuronal tissues [RGD, Feb 2006]



[View online »](#)