

Product datasheet for **RR203125L4V**

Slc1a3 (NM_019225) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc1a3 (NM_019225) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc1a3
Synonyms:	EAAT1; GLAST; GLAST-1; GluT-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019225
ORF Size:	1629 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203125).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019225.1 , NP_062098.1
RefSeq Size:	3864 bp
RefSeq ORF:	1632 bp
Locus ID:	29483
UniProt ID:	P24942
Cytogenetics:	2q16
Gene Summary:	carries out Na ⁺ dependent transport of glutamate and aspartate; may regulate neurotransmitter concentration in the central nervous system [RGD, Feb 2006]



[View online »](#)