

## Product datasheet for **RR203138L3V**

### **S100a4 (NM\_012618) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	S100a4 (NM_012618) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	S100a4
Synonyms:	18A2; 42A; CAPL; MTS1; P9ka; PEL98; RNP9KA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012618
ORF Size:	303 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203138).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012618.2</a> , <a href="#">NP_036750.1</a>
RefSeq Size:	486 bp
RefSeq ORF:	306 bp
Locus ID:	24615
UniProt ID:	<a href="#">P05942</a>
Cytogenetics:	2q34
Gene Summary:	may play a role in regulation of cell growth and differentiation [RGD, Feb 2006]



[View online »](#)