

## Product datasheet for **RR203142L4V**

### Cttn (NM\_021868) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cttn (NM_021868) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cttn
Synonyms:	Cttnb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021868
ORF Size:	1527 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203142).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021868.2</a> , <a href="#">NP_068640.2</a>
RefSeq Size:	2951 bp
RefSeq ORF:	1530 bp
Locus ID:	60465
UniProt ID:	<a href="#">Q66HL2</a>
Cytogenetics:	1q41-q42
Gene Summary:	component of clathrin-coated pits; involved in receptor-mediated endocytosis [RGD, Feb 2006]



[View online »](#)