


## Product datasheet for RR203163L3

### Mob4 (NM\_133528) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Mob4 (NM_133528) Rat Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Mob4
Synonyms:	Mobkl1; Mobkl3; phocein; Phocn; Prei3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203163).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_133528
ORF Size:	675 bp



[View online >](#)

