

Product datasheet for **RR203168L3V**

Rab3d (NM_080580) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rab3d (NM_080580) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rab3d
Synonyms:	MGC93243
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_080580
ORF Size:	657 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203168).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080580.2 , NP_542147.1
RefSeq Size:	2226 bp
RefSeq ORF:	660 bp
Locus ID:	140665
UniProt ID:	Q63942
Cytogenetics:	8q13
Gene Summary:	may regulate membrane fusion in Fc epsilonRI-stimulated mast cell exocytosis [RGD, Feb 2006]



[View online »](#)