

Product datasheet for **RR203174**

Nat2 (NM_053854) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nat2 (NM_053854) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nat2
Synonyms: Nat2a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR203174 representing NM_053854
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGACATTGAAGCATACTTTGAAAGAATTGGTTATCAGAGTTCAGGAACAACTGGACTTGAAGAAT
TAACTGAAATCTTCAACACCAGATACGGGCTATCCCTTTGAAAACCTGAACATCCATTGCGGGGAATC
CATGGAGCTGAATTTAGAAGTCATCTTTGATCAAGTTGTGAGGAAGAAGCGTGGTGGATGGTGTCTCCAG
GTTAATCATCTGCTGACTGGGCTCTGACCAAAATGGGCTTTGAAGCCACAATGTTGGGAGGGTATGTCT
TTAACTCCAGCCAATAAGTACAGCAGTGGCATGATTCACTTACTAGTGCAGGTGACCCCTCAGTGGAAA
GGACTACATTGTTGATGCTGGGTTTGGACGTTCCCTACCAGATGTGGGAGCCTCTGGAATTAACATCAGGG
AAGGATCAGCCTCAGGTGCCTGCCATCTTTCTGTTGACAGAAGAGAATGGAACCTGGTACTTGGACCAA
TCAGAAGAGAGCAGTATGTTCCAAACCAAGAGTTTGTAACTCAGACCTCCTTAAAAGAACAAATATCG
AAAAATCTATTCCTTTACTCTTGAGCCCCGCACTATTGAGGATTTGAGTCCATTAATACCTATCTTCAG
ACATCGCCAGCATCCTTGTACAAGTAAATCATTTTGTCTTGCAGACCCTAGAAGGGGTTCACTGTT
TGGTTGGCTCCACCCTCACTTATAGGAGATTCAGTTATAAGGACAATATAGATCTTGTAGAGTTTAAAGAG
TCTGACTGAGGAAGAAATAGAAGACGTAAGTAACTATATTTGGGGTTCTTTAGAGAGAAAGCTTGTG
CCTAAACATGGTATCGATTTTTTACCATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR203174 representing NM_053854
 Red=Cloning site Green=Tags(s)

MDIEAYFERIGYQSSRNKLDLEELTEILQHQIRAI PFENLN IHCGESMELNLEVIFDQVVRKKRGGWCLQ
 VNHLLYWALTKMGFEATMLGGYVFNT PANKYSSGMIHLLVQVTL SGKDYIVDAGFGRSYQMWEPLEL TSG
 KDQPQVPAIFRL TEENGTWYLDQIRREQYV PNQEFVNSDLLEKNKYRKIYSFTLEPRTIEDFESINTYLQ
 TSPASLFTSKSFCSLQTLEGVHCLVGSTLYRRFSYKDNIDLVEFKSLTEEEIEDVLKTI FGVSLERKLV
 PKHGDRFFTI


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_053854

ORF Size: 870 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_053854.1 , NP_446306.1
RefSeq Size:	1385 bp
RefSeq ORF:	873 bp
Locus ID:	116632
UniProt ID:	P50298
Cytogenetics:	16p14
MW:	33.8 kDa
Gene Summary:	enzyme that acetylates only aryl amines; may have a role in carcinogenesis [RGD, Feb 2006]