

## Product datasheet for **RR203183**

### Camp (NM\_001100724) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Camp (NM\_001100724) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Camp  
**Synonyms:** CRAMP  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RR203183 representing NM\_001100724  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCAACCCCATAGGGACGTCCTTCCCTGTGGCGCTCACTGTCACTGCTATTGCTCCTGGGCCTTGGGT  
 TGCCTCTAGCCGTTTCCAGACCCTCAGCTACAGGAGGCTGTACTCCGTGCTGTGGATGACTTCAACCA  
 GCAGTCTTTGGACCAATCTCTACCGTCTCCTGGACCTGGATTCTGAGCCCCAAGGGGATGAGGACCA  
 GATACTCCCAAGTATGTGAGGTTCCGAGTGAAGGAGACTGTCTGTAGCAAGGCATCACAGCAGTTACCTG  
 AGCAATGTGCCTTCAAGGAACAGGGGTAGTGAAGCAGTGTATGGGGACAGTCACCCTGAACCAGGCTGC  
 AGAGTCTTTTGACATCAGTTGTGATGCGCCTGGTATACAGCCTTTCCGGTTCAAGAAAATTTCCCGGCTG  
 GCTGGACTCCTCCGAAAAGGTGGAGAGAAGTTTGGTAAAAGCTTAGGAAAATCGGCCAGAAAATTAAGG  
 ATTTTTCCAGAACTTGACACTGAAATAGAGCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR203183 representing NM\_001100724  
 Red=Cloning site Green=Tags(s)

MQPHRDVPSLWRSLSLLLLLGLGLPLAVSQTLSYREAVLRAVDDFNQQSLDTNLYRLDLLDSEPQGDDEP  
 DTPKYVRFVKETVCSKASQQLPEQCAFKEQGVVKQCMGTVTLNQAAESFDISCDAPGIQPFRRFKISRL  
 AGLLRKGGEKFGEKLRKIGQKIKDFFQKLAPEIEQ

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

## Cloning Scheme:



ACCN: NM\_001100724

ORF Size: 525 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.  
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.  
3. Close the tube and incubate for 10 minutes at room temperature.  
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.  
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001100724.1](#), [NP\\_001094194.1](#)

RefSeq Size: 586 bp  
 RefSeq ORF: 528 bp  
 Locus ID: 316010  
 Cytogenetics: 8q32  
 MW: 19.8 kDa

## Product images:



Circular map for RR203183