

Product datasheet for **RR203200L3V**

Arc (NM_019361) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arc (NM_019361) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Arc
Synonyms:	rg3.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019361
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203200).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019361.1 , NP_062234.1
RefSeq Size:	3032 bp
RefSeq ORF:	1191 bp
Locus ID:	54323
UniProt ID:	Q63053
Cytogenetics:	7q34
Gene Summary:	may play a role in neuronal plasticity, long-term potentiation, and memory consolidation [RGD, Feb 2006]



[View online »](#)