

Product datasheet for **RR203204L3V**

Hdgfrp3 (Hdgfl3) (NM_145785) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hdgfrp3 (Hdgfl3) (NM_145785) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hdgfl3
Synonyms:	Hdgfrp3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145785
ORF Size:	606 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203204).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145785.1 , NP_665728.1
RefSeq Size:	2750 bp
RefSeq ORF:	609 bp
Locus ID:	252941
UniProt ID:	Q923W4
Cytogenetics:	1q31
Gene Summary:	Enhances DNA synthesis and may play a role in cell proliferation.[UniProtKB/Swiss-Prot Function]



[View online »](#)