

Product datasheet for **RR203222L3V**

P2rx2 (NM_053656) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	P2rx2 (NM_053656) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	P2rx2
Synonyms:	P2X2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053656
ORF Size:	1416 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203222).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053656.3
RefSeq Size:	1833 bp
RefSeq ORF:	1419 bp
Locus ID:	114115
UniProt ID:	P49653
Cytogenetics:	12q16
Gene Summary:	subunits assemble into ATP-gated ion channels involved in neuritogenesis and cell survival [RGD, Feb 2006]



[View online »](#)