

Product datasheet for **RR203229**

Nfic (NM_031567) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nfic (NM_031567) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nfic
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR203229 representing NM_031567
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTATTCCTCCCCGCTCTGCCTGACCCAGGATGAGTTCACCCATTCATCGAGGCTCTGCTGCCGCATG
TCCGTGCCTTCGCTACACCTGGTTCAACCTGCAGGCGCAAGCGCAAGTACTTCAAGAAGCACGAGAA
GCGCATGTGCAAGGACGAGGAGCGCGGTGAAGGATGAACTGCTGGGTGAAAAGGCCGAGGTGAAGCAG
AAGTGGGCGTCACGGCTGCTGGCCAAGCTGCGCAAGGACATCCGGCCGAGTGCCCGAAGACTTCGTAC
TGGCTGTGACAGGCAAGAAGGCACCAGGCTGTGTGCTGTCCAACCCGGACCAGAAGGGCAAGATGCGGCG
CATCGACTGCCTGCGCCAGGCTGACAAAGTATGGCGCCTGGACCTGGTTCATGGTTCATCCTGTTCAAGGGT
ATCCCACTGGAGAGCACGGACGGCGAGCGGCTGGTGAAGGCCGCTGCCTGCGCGCACCCCGTGTGTGTG
TGCAGCCGACCATATTGGCGTGGCCGTCAAGGAGCTGGACCTGTACCTGGCCTACTTCGTGCGGGAGCG
CGATGCGGAGCAGAGCAGCAGTCCCGGACAGGAGTGGGCTCAGACCAGGAGGACAGCAAGCCATTACA
CTGGACACGACAGACTTCCAGGAGACTTCGTACATCTGGCGTGTTCAGTGTACGGAGCTGATCCAGG
TGTCTCGGACACCTGTGGTGACCGGAACCGGGCCAACTTCTCCCTGGGAGAGCTGCAGGGCCACCTGGC
ATACGACCTGAACCCAGCCAGCGCTGGCATGAGAAGGACTCTACCCAGTACTTCCCTCCAGTGGGAGCAAG
CGGCACAAGTCCGGCTCGATGGAGGAGGACGTGGACACGAGCCCTGGCGGAGACTACTACACATACCCCA
ACTCGCCACGAGTAGCAGCCGCAACTGGACGGAAGACATGGAAGGAGGCATCTCATCCCGGTGAAGAA
GACAGAGATGGACAAATCTCAATTAACAGCCATCCCCCAGGACTCCCCCGGCTCTCCAGCTTACACA
CAGCACACAGGCCCGTCATCGCGGTACACAGCGGCATCGCGGGAGCCCCACCCGTCCTCAGCCCTGC
ACTTCCCTGCCACGCCATCCTGCCGACAGAGCTCCACCTACTTCCCCCACTGCCATCCGGTACCC
GCCACATCTCAACCCAGGACCACTCAAGGATCTGGTGTCCCTGGCCTGCGACCCGGCCACCCAGCAG
CCAGGACCGCTGCTCTCCGCCGACGTCCCCTGCAACAGTCCCTTTGTGGAT

**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA**




Protein Sequence: >RR203229 representing NM_031567
 Red=Cloning site Green=Tags(s)

MYSSPLCLTQDEFHPFIEALLPHVRAFAYTWFNLQARKRKYFKKHEKRMKDEERAVKDELLGEKAEVKQ
 KWARLLAKLRKDIRPECREDFLAVTGKKAPGCVLSNPDQKGMRRIDCLRQADKVVRLDLVMVILFKG
 IPLESTDGERLVKAAACAHPVLCVQPHHIGVAVKELDLAYFYFRERDAEQSSSPRTGVGSDQEDSKPIT
 LDTTDFQESFVTSQVFSVTELIQVSRTPVVTGTGPNFSLGELQGHLAYDLNPASAGMRRTLPSTSSSGSK
 RHKSGSMEEDVDTSPGGDYTSPNSPTSSRNWTEDEMEGGISSPVKKTMDKSPFNPSPPQDSPRLSSFT
 QHHRPVIIVHSGIARSPHPSSALHFPATPILPQTASTYFPHTAIRYPPHLNPQDPLKDLVSLACDPATQQ
 PGPPALRPARPLQTVPLWD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_031567
ORF Size:	1317 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_031567.1 , NP_113755.1
RefSeq Size:	1565 bp
RefSeq ORF:	1320 bp
Locus ID:	29228
Cytogenetics:	7q11
MW:	48.7 kDa
Gene Summary:	may regulate gene expression in olfactory neurons [RGD, Feb 2006]