

Product datasheet for RR203242L4

Cckar (NM_012688) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cckar (NM_012688) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Cckar
Synonyms:	Cck-ar
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203242).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 24889
UniProt ID: [P30551](#)
Cytogenetics: 14q11
Gene Summary: G-protein coupled receptor activated by cholecystokinin; may be involved in food intake and insulin regulation [RGD, Feb 2006]